

Cortical Thickness

Pipeline v1.0 Overview

Usage

Setup

SET BASE_PATH indicating the path to CorticalThickness folder.

Example:

```
BASE_PATH=${1:-"/neuro/labs/grantlab/research/MRI_processing/jose.cisneros/CorticalThickness"}
```

Usage without Docker.

```
(base) [jose.cisneros@hanyang:x86_64-Linux]...rticalThickness$>python code/corticalThickness.py -h
usage:  ===== Cortical Thickness Pipeline by Jose Cisneros (May 5, 2022 ver.1) =====

    [-h] -ca CASE [-bp BASE_PATH] [-im IN_MRI] [-is IN_MRI_SEG] [-se ENABLE_SURFACE_EXTRACTION]
    [-ir ENABLE_INTENSITY_REFINEMENT] [-do OUTSIDE_DOCKER] [-rp RESULTS_PREFIX] [-icm INTENSITY_CLUSTERING_METHOD]
```

optional arguments:

```
-h, --help            show this help message and exit
-ca CASE, --case CASE
                        Subject id
-bp BASE_PATH, --base-path BASE_PATH
                        Path containing code & resources.
-im IN_MRI, --input-mri IN_MRI
                        Input MRI .nii file path.
-is IN_MRI_SEG, --input-segmented IN_MRI_SEG
                        Input MRI segmented .nii file path.
-se ENABLE_SURFACE_EXTRACTION, --surface-extraction ENABLE_SURFACE_EXTRACTION
                        Enable Surface Extraction (default: True)
-ir ENABLE_INTENSITY_REFINEMENT, --intensity-refinement ENABLE_INTENSITY_REFINEMENT
                        Enable Intensity Clustering for CP external boundary refinement. (default: True)
-do OUTSIDE_DOCKER, --outside-docker OUTSIDE_DOCKER
                        Flag indicating if script running outside docker. (default: True)
-rp RESULTS_PREFIX, --results-prefix RESULTS_PREFIX
                        Prefix for Results folder. (default: )
-icm INTENSITY_CLUSTERING_METHOD, --intensity-clustering-method INTENSITY_CLUSTERING_METHOD
                        Soft Clustering Method use for Intensity Refinement. Options: GMM, FCM, sFCM (default: sFCM)
```

Output Folder: {BASE_PATH}/Results/{CASE}

Example:

```
${BASE_PATH}/code/corticalThickness.py \
-ca FCB028 \
-im ${BASE_PATH}/Samples/FCB028/recon_to31.nii \
-is ${BASE_PATH}/Samples/FCB028/segmentation_to31_final.nii
```

Usage with Docker.

```
(base) [jose.cisneros@hanyang:x86_64-Linux]...rticalThickness$>python code/corticalThicknessDocker.py -h
usage:  ===== Cortical Thickness Pipeline using Docker by Jose Cisneros (May 9, 2022 ver.1) =====

    [-h] -ca CASE [-bp BASE_PATH] -im IN_MRI -is IN_MRI_SEG -o OUTPUT [-se ENABLE_SURFACE_EXTRACTION]
    [-ir ENABLE_INTENSITY_REFINEMENT] [-icm INTENSITY_CLUSTERING_METHOD]

optional arguments:
  -h, --help                show this help message and exit
  -ca CASE, --case CASE      Subject id
  -bp BASE_PATH, --base-path BASE_PATH
                              Path containing code & resources.
  -im IN_MRI, --input-mri IN_MRI
                              Input MRI .nii file with path.
  -is IN_MRI_SEG, --input-segmented IN_MRI_SEG
                              Input MRI segmented .nii file with path.
  -o OUTPUT, --output-folder OUTPUT
                              Output folder containing all generated files.
  -se ENABLE_SURFACE_EXTRACTION, --surface-extraction ENABLE_SURFACE_EXTRACTION
                              Enable Surface Extraction
  -ir ENABLE_INTENSITY_REFINEMENT, --intensity-refinement ENABLE_INTENSITY_REFINEMENT
                              Enable Intensity Clustering for CP external boundary refinement.
  -icm INTENSITY_CLUSTERING_METHOD, --intensity-clustering-method INTENSITY_CLUSTERING_METHOD
                              Soft Clustering Method use for Intensity Refinement. Options: GMM, FCM, sFCM
```

Example:

```
$(BASE_PATH)/code/corticalThicknessDocker.py \
  -ca FCB028 \
  -im $(BASE_PATH)/Samples/FCB028/recon_to31.nii \
  -is $(BASE_PATH)/Samples/FCB028/segmentation_to31_final.nii \
  -o $(BASE_PATH)/Results
```

Useful information about Docker.

Setup:

- Option1: Get Docker Image building it with Dockerfile
 - `docker build -t cortical-thickness ${BASE_PATH} -f ${BASE_PATH}/deploy/Dockerfile`
- Option2: Get latest Docker Image published.
 - `docker pull ghcr.io/josecisneros001/fnndsc-cortical-thickness`
 - `docker tag ghcr.io/josecisneros001/fnndsc-cortical-thickness cortical-thickness`

Publish new Docker Image


- `docker login ghcr.io`
- # Build your image if it doesn't exist yet.
 - `docker build -t ghcr.io/josecisneros001/FNNDSC-Cortical-Thickness .`
- # Tag it if already exists.
 - `docker tag cortical-thickness ghcr.io/josecisneros001/fnndsc-cortical-thickness`
- `docker push ghcr.io/josecisneros001/fnndsc-cortical-thickness`

Code - Pipeline


Processing Step

Segment and Join Cerebral Exterior Labels 1 & 42. Obtained from deep learning segmentation.


`${TARGET_DIR}/${CASE}/temp/cerebral_ext.mnc`





```
# Dilation Cerebral Exterior
${TARGET_DIR}/${CASE}/temp/cerebral_ext_d.mnc
```




Segment and Join Cerebral Interior Labels 160 & 161. Obtained from deep learning segmentation.
\${TARGET_DIR}/\${CASE}/temp/cerebral_int.mnc




Dilation Cerebral Interior
\${TARGET_DIR}/\${CASE}/temp/cerebral_int_d.mnc



Binarize all segmentations [from 1 to 161] Obtained from deep learning segmentation.
\${TARGET_DIR}/\${CASE}/temp/initial_segmentations.mnc




```
# Dilation Initial Segmentations
${TARGET_DIR}/${CASE}/temp/initial_segmentations_d.mnc
```




```
# Remove Cerebellum for next steps.  
${TARGET_DIR}/${CASE}/input/input_mri_processed.mnc
```

Before:




After:



Code - Pipeline

Skeletonization Step

```
# GM External Boundary - 1 voxel apart.  
${TARGET_DIR}/${CASE}/temp/cerebral_ext.mnc -> White (Input)  
${TARGET_DIR}/${CASE}/temp/gm_ext.mnc -> Blue (Output)
```



Join GM and WM, labeled, not binarized


\${TARGET_DIR}/\${CASE}/temp/gray.mnc (Input)

\${TARGET_DIR}/\${CASE}/temp/cerebral_int.mnc (Input)

\${TARGET_DIR}/\${CASE}/temp/wm_and_gm.mnc -> White (GM) & Pink (WM) (Output)

GM with GM/CSF outer boundary + WM



Important detail to consider for the skeleton input. GM and WM need to have different labels, that's why the output has two different colors (Pink & White).



```
#### Prepare Skeleton Gray Scale input.
#### 11-Outside GM, 0 - WM, 255 - GM.
${TARGET_DIR}/${CASE}/temp/cerebral_ext.mnc (Input)
${TARGET_DIR}/${CASE}/temp/cerebral_int.mnc (Input)
${TARGET_DIR}/${CASE}/temp/gm_grayscale.mnc -> Blue (GM), Dark (WM) & lightBlue (outside) (Output)
```



CSF obtained with CP dilation + GM

Important detail that lightblue is everything the skeleton don't need to consider. First Image doesn't hide lightblue to show that there are parts of the brain not considered in skeleton.



Skeleton

```
${TARGET_DIR}/${CASE}/temp/gm_grayscale.mnc (Image 2 Blue) (Input)  
${TARGET_DIR}/${CASE}/temp/wm_and_gm.nii (Image 1 White & Pink) (Input)  
${TARGET_DIR}/${CASE}/temp/skeleton_1_corr.nii -> (Image 1&2 Red) (Output)
```





Join Skeleton with GM external boundary.Skeleton follows External boundary limited by DeepLearning in a 100%.

\${TARGET_DIR}/\${CASE}/temp/gm_ext.mnc (Input First Image)

\${TARGET_DIR}/\${CASE}/temp/skeleton_1_corr.nii (Input Second Image Blue)

\${TARGET_DIR}/\${CASE}/temp/skeleton_corr.nii -> (Output Second Image Blue + Red)



Code - Pipeline


Intensity Refinement (Optional)

Clustering Input Volume: GM -> 1 extra outer voxel and without 1 inner voxel.

\${TARGET_DIR}/\${CASE}/temp/clustering_input.mnc (Blue)

Purple - GM.

Input for Soft clustering method,
divide voxels in two groups (Bright
& Dark).




Current Pipeline has three soft clustering methods implemented. Choose by cmdline argument.

- Gaussian Mixture Model (GMM)
- FCM (Fuzzy c-Means)
- sFCM (Spatial Fuzzy c-Means)


GMM & FCM are based only in Intensity, sFCM considers neighboring intensity, for this reason it has a better performance against noise and tends to have all voxels connected.

sFCM example:




Join Clustering Output with Skeleton limiting inner boundary to skeleton.


\${TARGET_DIR}/\${CASE}/temp/ps2_csf.mnc (Red)





```
#### Get inner boundary.  
${TARGET_DIR}/${CASE}/temp/ps2_csf.mnc (Red)  
Red Inner Boundary:  
${TARGET_DIR}/${CASE}/output/ps2.nii (Blue)  
also called as:  
${TARGET_DIR}/${CASE}/output/skeleton_output.mnc
```



Skeleton



Skeleton + Intensity Refinement




Code - Pipeline

Surface Extraction


White Inner Extraction

```
${TARGET_DIR}/${CASE}/surfaces/lh.smoothwm.native.obj  
${TARGET_DIR}/${CASE}/surfaces/lh.smoothwm.native.asc
```


```
${TARGET_DIR}/${CASE}/surfaces/rh.smoothwm.native.obj  
${TARGET_DIR}/${CASE}/surfaces/rh.smoothwm.native.asc
```



```
#### Generate Laplace Field for outer extraction.  
${TARGET_DIR}/${CASE}/surfaces/laplacian_to31_left_fixed.mnc
```



```
#### Generate Laplace Field for outer extraction.  
${TARGET_DIR}/${CASE}/surfaces/laplacian_to31_right_fixed.mnc
```




Expand White Matter till Laplace Field, CP Surface Extraction

\${TARGET_DIR}/\${CASE}/surfaces/lh.pial.native.obj

\${TARGET_DIR}/\${CASE}/surfaces/lh.pial.native.asc

\${TARGET_DIR}/\${CASE}/surfaces/rh.pial.native.obj

\${TARGET_DIR}/\${CASE}/surfaces/rh.pial.native.asc;



Code - Pipeline

Morphometrics Measurement

Cortical Thickness Calculations with three methods (same as Mengting Pipeline) (3 .txt files per side generated.)

- tlink use t-link method for thickness
- tnear use t-near method for thickness
- tlaplace use t-Laplace method for thickness

Depth Potential Calculations. (1 .txt file per side generated.)

Results/FCB028/morphometrics/*

▼ morphometrics
≡ native_rms_tlaplace_10mm_left.txt
≡ native_rms_tlaplace_10mm_right.txt
≡ native_rms_tlink_10mm_left.txt
≡ native_rms_tlink_10mm_right.txt
≡ native_rms_tnear_10mm_left.txt
≡ native_rms_tnear_10mm_right.txt
≡ wm_sulcaldepth_left.txt
≡ wm_sulcaldepth_right.txt

Format:

1	15.021
2	27.0097
3	28.6648
4	27.4705
5	47.6426
6	55.6526
7	46.3256
8	54.2754
9	28.0118
10	33.8271
11	32.9404
12	24.1953
13	22.9465
14	18.7177
15	10.2492
16	52.3813
17	39.4696
18	33.0017
19	36.0228
20	40.8147
21	18.3901
22	56.9144
23	36.5423
24	48.1512
25	26.8554
26	40.9527
27	25.3158
28	44.0519
29	39.5022
30	36.461
31	36.4782
32	34.0717
33	17.2389
34	23.7858
35	7.52124
36	43.3553
37	61.2423
38	46.7108
39	47.5814
40	32.8099